



FEVERSCAN Disposable Forehead Thermometer

The FEVERSCAN Disposable Forehead Temperature Monitor for Continuous Digital Temperature Event Monitoring For The Perioperative Environment.


These liquid crystal thermometers are to be used on the forehead area only to estimate the core body temperature. This is an ideal way to check for temperature changes. These strips are meant to be used as an unbreakable, quick reading device to be used between clinical thermometer readings and to monitor trends.

This series is available in two models:

- FSP430 94°F to 104°F
- FSP431 35°C to 40°C

DESCRIPTION

- DIM 1.8" W x 0.4"H
- Nominal 0.01" thickness
- 6 increments in °F or °C



100 per box

TO READ

- Remove backing. Apply to center of forehead. Skin surface should be clean and dry
- Read temperature after colors stop changing. Typically 15 seconds
- GREEN. Correct reading. If GREEN not visible, subtract 1°F from highest lit channel to obtain correct reading.

PACKAGING

- 100 ea in Lot numbered box

COMPONENTS

- Hypoallergenic pressure sensitive acrylate on 3 mil polyethylene film
- Polyester substrate
- Non-toxic mixture of liquid crystal esters

PRODUCT PERFORMANCE

- Accuracy $\pm 1.0^{\circ}\text{F}$ / 0.6°C
- Use to monitor temperature trends between clinical thermometer readings
- Core adjusted with 4.0°F / 2.2°C offset

SHELF LIFE: Guaranteed 1 years from date of shipment if maintained in a cool (less than 21°C / 70°F), and dry less than 50% RH) environment.

SHIPMENT: Irreversible temperature recording labels with ascending response temperatures of 130°F (54°C) or lower must ship via 2nd Day Air or faster and may seasonally incur an additional charge for insulated packaging.

WARRANTY: The following is made in lieu of all warranties, expressed or implied, including the implied warranties of merchantability and fitness of purpose: Seller's and manufacturers only obligation shall be to replace such quantity of the product proved by Telatemp Corporation to be defective.

LIMITATION OF LIABILITY: All physical properties, statements, and recommendations are either based on tests we believe to be reliable or our experience, but they are not guaranteed. Telatemp Corporation recommends each user determine the suitability of the products for the intended application. Telatemp Corporation shall not be liable for direct, indirect, special, incidental or consequential damages based upon breach of warranty, breach of contract, negligence, strict liability or any other legal theory.

**Prices and specifications are subject to change without notice.*